

NEW

S-CAPE®

## EXTERNALDISPLAY

24" | 27" | 32"



typical configuration

Ready for the **CX-PRIME** integration environment.

## Pendant displays for convenient image viewing.

S-CAPE® made it their business to adapt the OR-Situation to the individual requirements and special tasks of each OR-Team in the context of high-tech image viewing.

### Fits your needs with:

- Fanless, enclosed display with ESG safety glass
- IP 65 protection class
- Medical power supply compliant to the required device leakage currents
- Compliant to hygienic requirements for the operation in sterile environments
- Utilizes the performance of the S-CAPE® **MULTICONSOLE**
- Digital image transmission in HD quality

Further information about the **EXTERNALDISPLAY** and the **CX-PRIME** integration solution:

S-CAPE® GmbH  
Komturstraße 18a, 12099 Berlin  
T: +49 30 7130297-0  
E: info@s-cape.com  
www.s-cape.com

### technical data

#### housing:

- + aluminium housing: **powder-coated RAL 9010**

#### EXTERNALDISPLAY 24"

- + Technology: **LED**
- + Contrast Ratio: **3000 : 1**
- + Brightness: **250 cd/m<sup>2</sup>**
- + Response Time: **8.5 ms**
- + Vertical/Horizontal Viewing Angle: **178/178**
- + Input Modules (one of the below):  
**3G/HD-SDI or DVI Single-Link or CAT-5e / CAT-6**

#### EXTERNALDISPLAY 27"

- + Technology: **LED**
- + Contrast Ratio: **3000 : 1**
- + Brightness: **350 cd/m<sup>2</sup>**
- + Response Time: **7 ms**
- + Vertical/Horizontal Viewing Angle: **178/178**
- + Input Modules (one of the below):  
**3G/HD-SDI or DVI Single-Link or CAT-5e / CAT-6**

#### EXTERNALDISPLAY 32"

- + Technology: **LED**
- + Contrast Ratio: **1600 : 1**
- + Brightness: **450 cd/m<sup>2</sup>**
- + Response Time: **8 ms**
- + Vertical/Horizontal Viewing Angle: **178/178**
- + Input Modules (one of the below):  
**3G/HD-SDI or DVI Single-Link or CAT-5e / CAT-6**