

NEW

S-CAPE 

TELEPORT



Ready for **CX-PRIME SHARE**

with wired or wireless network access

Enabling **remote** collaboration **in real-time** surgical setting.

TELEPORT's mobile, cleanroom-compliant design can establish remote presence at the point of care within minutes.

Tailored to your needs with:

- Interoperable with endpoints of all major telepresence vendors
- Compatible with specialized telemedicine solutions from **CX-PRIME SHARE**
- Sealed cable management and monitor enclosure for IP65-compliant infection control in mission-critical settings
- Optional display sizes and technologies
- Optional add-ons include battery power, wireless network access

Further information about the **TELEPORT** and the **CX-PRIME** integration solution:

S-CAPE® GmbH
Komturstraße 18a, 12099 Berlin
T: +49 30 7130297-0
E: info@s-cape.com
www.s-cape.com

technical data

monitor options

- + 24"-LED-Monitor
Resolution: 1920x1080 (FullHD)
Brightness: 250 cd/m²
Contrast: 4000 : 1
Optional: Touch Upgrade
- + 27"-LED-Monitor
Resolution: 1920x1080 (FullHD)
Brightness: 350 cd/m²
Contrast: 3000 : 1

powered by:

- + Mainboard with INTEL® Core™ i5 processor
- + 4 GB RAM (expandable up to 32 GB)
- + LAN RJ 45 Gigabit (10/100/1000 mb/s)
- + 500 GB SATA HDD or other options possible

features:

- + Supports single or dual monitor configuration with IP65 adherence
- + Advanced mobility enables usage across different procedure rooms
- + Fully enclosed system offers easier cleaning, reduces infection control concerns and prevents ingress of liquids
- + Antimicrobial and easy-clean finish
- + Lockable casters